



Recombinant Chicken Protein cereblon (CRBN)

Product Code	CSB-EP005944CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0CF65
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	MAAEEGGDGR RNMGNPPPPA PAESEEEEDDN EMEVEDQDGK EAEKPNMINF DTS�PTSHMY LGSDMEEFHG RTLHDDDDSCQ VIPVLPHEMV MLIPGQTLPL QLFHPQEVSM VRNLIQKDRT FAVLAYSNVR EREAHFHTTA EIIAYREEQE YGIETVKVKA IGRQRFKVL E IRTQSDGIQQ AKVQILPERV LPSTMSAVQL QSLSRRHIFP SSKPKVWQDR AFRQWWQKYQ KRKFHCASLT SWPPWLYSLY DAETLMERVK RQLHEWDENL KDESŁPTNPI DFSYRVAACL PIDDALRIQL LKIGSAIQRL RCELDIMNKC TSLCCKQCQD TEITTKNEIF SLSLCGPMAA YVNPHGYIHE TLTVYKACNL NLSGRPSTEH SWFPGYAWTI AQCRICGNHM GWKFTATKGD MSPQKFWGLT RSALLPRIPE AEDELGHDRS PLLCL
Source	E.coli
Target Names	CRBN
Protein Names	Recommended name: Protein cereblon
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.