



Recombinant Chicken Protein C-ets-2 (ETS2)

Product Code	CSB-MP007853CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P10157
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MSEFAIRNMD QVAPVSNMYR GMLKRQPAFD TFDSSNSLFA GYFLSLNEDQ TLQEVPTGFD STSYESNCE LPLLTPCSKA VMSQALKDTF SGFTKEQCRL GIPNNPWLWT EQHVCQWLAW ATNEFSLANV NIHQFLMSGQ DLCNLGKERF LELAPDYVGD ILWEHLEQMI KDSQEKTQDQ YVESSHLTSV PHWVNNNSLT VNVDQTPYGI QMPGYPKALS YPKPNLLSDI CQTSTGPNLL SPEQDFSLFP KTQVDAVSVN YCTVNQDFTR SNLNLIDNS GKLREHESSE SGAESYESSD SMLQSWNSQS SLVDLQRVPS YESFEDDCSQ SLCMSKPTMS FKDYIQDRSD PVEQGKPVIP AAILAGFTGS GPIQLWQFLL ELLTDKSCQS FISWTGDGWE FKLADPDEVA RRWGRRKNKP KMNYEKLSRG LRYYYDKNII HKTSGKRYVY RFVCDLQNLG YTAEEELHAM LGVQPDTED
Source	Mammalian cell
Target Names	ETS2
Protein Names	Recommended name: Protein C-ets-2
Expression Region	1-479
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.