



Recombinant Gossypium barbadense NAD (P)H-quinone oxidoreductase subunit K, chloroplastic

Product Code	CSB-BP371605GEP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0ZZ39
Product Type	Recombinant Protein
Immunogen Species	Gossypium barbadense (Sea-island cotton) (Egyptian cotton)
Purity	≥85% (SDS-PAGE)
Sequence	MNSIEFPLLD RTTQNSVIST TLNDLSNWSR LSSLWPLLYG TSCCFIEFAS LIGSRFDFDR YGLVPRSSPR QADLILTAGT VTMKMAPSLV RLYEQMPEPK YVIAMGACTI TGGMFSTDSY STVRGVDKLI PVDVYLPGCP PKPEAVIDAI TKLRKKISRE IYEDRIRSQQ GDRCFTTNHK FCLVRSTRTG NYNQGLLYQP PSTSEIPPET FFNYKGSLSLSS HELVN
Source	Baculovirus
Target Names	ndhK
Protein Names	Recommended name: NAD(P)H-quinone oxidoreductase subunit K, chloroplastic EC= 1.6.5.- Alternative name(s): NAD(P)H dehydrogenase subunit K NADH-plastoquinone oxidoreductase subunit K
Expression Region	1-225
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.