



Recombinant *Gossypium barbadense* 30S ribosomal protein S3, chloroplastic (rps3)

Product Code	CSB-BP371616GEP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0ZZ73
Product Type	Recombinant Protein
Immunogen Species	<i>Gossypium barbadense</i> (Sea-island cotton) (Egyptian cotton)
Purity	>85% (SDS-PAGE)
Sequence	MGQKINPLGF RLGTTQSHHS LWFAQPKKYS EGLQEDKKIR DCIKNYVQKN TRLSSGVEGI ARIEIQRKLD LIQVIIMGF PKLLIEDKPR KLEELQMNVQ KELNCMNRKL NIAITRIGNP YGHPNILAEF IAGQLKNRVS FRKAMKKAIE LTEQADTKGI QIQIAGRIDG KEIARVEWIR EGRVPLQTIG AKIEYCSYRV RTIYGVLGIK IWIFIDEE
Source	Baculovirus
Target Names	rps3
Protein Names	Recommended name: 30S ribosomal protein S3, chloroplastic
Expression Region	1-218
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.