



Recombinant Allantoinase (allB)

Product Code	CSB-MP371730EJE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A1A8H7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O1:K1 / APEC
Purity	≥85% (SDS-PAGE)
Sequence	MSFDLIKNG TVILENEARV VDIKVDGKI AAIGQDLGDA KDVMASGLV VSPGMVDAHT HISEPGRSHW EGYETGTRAA AKGGITTMIE MPLNQLPATV DRASIELKFD AAKGKLTIDA AQLGGLVSYN IDRLHELDEV GVVGFKCFVA TCGDRGIDND FRDVNDWQFF KGAQKLGLG QPVLVHCENA LICDALGEEA KREGRVTAHD YVASRPVFTE VEAIRRVLYL AKVAGCRLHV CHVSSPEGVE EVTRARQEGQ DVTCECPHY FVLDTDQFEE IGTAKCSPP IRDLNQKGM WEKLFNGEID CLVSDHSPCP PEMKAGNIMK AWGGIAGLQS CMDVMFDEAV QKRGMSLPMF GKLMATNAAD IFGLQQKGRI APGKDAFVF IQPNSSYVLT NDDLEYRHKV SPYVGRTIGA RITKTILRGD VIYDIEQGFP VAPKGQFILK HQQ
Source	Mammalian cell
Target Names	allB
Protein Names	Recommended name: Allantoinase EC= 3.5.2.5 Alternative name(s): Allantoin-utilizing enzyme
Expression Region	1-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.