



Recombinant Pimelyl-[acyl-carrier protein] methyl ester esterase (bioH)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-EP371932EJE |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | A1AGT6 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O1:K1 / APEC |
| Purity | >85% (SDS-PAGE) |
| Sequence | MNNIWWQTKG QGNVHLVLLH GWGLNAEVWR CIDEELSSHF TLHLVDLPGF GRSRGFGALS LADMAEAVLR QAPDKAIWLG WSLGGLVASQ IALTHPERVQ ALVTVASSPC FSARDEWPGI KPDVLAGFQQ QLSDDFQRTV ERFLALQTMG TETARQDARA LKKTVLALPM PEVDVLNGGL EILKTVDLRL PLQNVSM PFL RLYGYLDGLV PRKVV PMLDK LWP HSESYIF AKA AHAPFIS HPVEFHLLLV ALKQRV |
| Source | E.coli |
| Target Names | bioH |
| Protein Names | Recommended name: Pimelyl-[acyl-carrier protein] methyl ester esterase EC= 3.1.1.85 Alternative name(s): Biotin synthesis protein BioH Carboxylesterase BioH |
| Expression Region | 1-256 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |