



Recombinant L-arabinose isomerase (araA)

Product Code	CSB-BP372148EJE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A1A7A9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O1:K1 / APEC
Purity	≥85% (SDS-PAGE)
Sequence	MTIFDNYEVW FVIGSQHLYG PETLRQVTQH AEHVVNALNT EAKLPCKLVL KPLGTTTPEI TAICRDANYD DRCAGLVVWL HTFSPAKMWI NGLTMLNKPL LQFHTQFNAA LPWDSIDMDF MNLNQTAHGG REFGFIGARM RQQHAVVTGH WQDKQAHERI GSWMRQAVSK QDTRHLKVCR FGDNMREVAV TDGDKVAAQI KFGFSVNTWA VGDLVQVVNS ISDGDVNALV DEYESCYTMT PATQIHGEKR QNVLEAARIE LGMKRFLEQG GFHAFTTTTFE DLHGLKQLPG LAVQRLMQQG YGFAGEGDWK TAALLRIMKV MSTGLQGGTS FMEDYTYHFE KGNDLVLGSH MLEVCPPIAV EEKPILDVQH LGIGGKDDPA RLIFNTQTGP AIVASLIDLG DRYRLLVNCI DTVKTPHSLP KLPVANALWK AQPDLPTASE AWILAGGAHH TVFSHALNLN DMRQFAEMHD IEITVIDNDT RLPAPKDALR WNEVYYGFRR
Source	Baculovirus
Target Names	araA
Protein Names	Recommended name: L-arabinose isomerase EC= 5.3.1.4
Expression Region	1-500
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.