



# Recombinant L-arabinose isomerase (araA)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-MP372148EJE  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | A1A7A9   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Escherichia coli O1:K1 / APEC  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)  |
| <b>Sequence</b>          | MTIFDNYEVW FVIGSQHLYG PETLRQVTQH AEHVVNALNT EAKLPCKLVL<br>KPLGTTTPEI TAICRDANYD DRCAGLVVWL HTFSPAKMWI NGLTMLNKPL<br>LQFHTQFNAA LPWDSIDMDF MNLNQT AHGG REFGFIGARM RQQHAVVTGH<br>WQDKQAHERI GSWMRQAVSK QDTRHLKVCR FGDNMREAV<br>TDGDKVAAQI KFGFSVNTWA VGDLVQVVNS ISDGDVNALV DEYESCYTMT<br>PATQIHGEKR QNVLEAARIE LGMKRFLEQG GFHAFTTTTFE DLHGLKQLPG<br>LAVQRLMQQG YGFAGEGDWK TAALLRIMKV MSTGLQGGTS FMEDYTYHFE<br>KGNDLVLGSH MLEVCPPIAV EEKPILDVQH LGIGGKDDPA RLIFNTQTGP<br>AIVASLIDLG DRYRLLVNCI DTVKTPHSLP KLPVANALWK AQPDLPTASE<br>AWILAGGAHH TVFSHALNLN DMRQFAEMHD IEITVIDNDT RLPAPKDALR<br>WNEVYYGFRR |
| <b>Source</b>            | Mammalian cell   |
| <b>Target Names</b>      | araA   |
| <b>Protein Names</b>     | Recommended name: L-arabinose isomerase EC= 5.3.1.4  |
| <b>Expression Region</b> | 1-500  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |