



# Recombinant Iron uptake system component EfeO (efeO)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP372214EJE   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | A1A9S2  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O1:K1 / APEC   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | ADIP QVKVTVTDKQ CEPMTITVNA GKTQFIIQNH SQKALEWEIL<br>KGVMVVEERE NIAPGFSQKM TANLQPGEYD MTCGLLTNPK GKLIVKGEAT<br>ADAAQSDALL SLGGAITAYK AYMMAETTQL VTDTKAFTDA IKAGDIEKAK<br>ALYAPTRQHY ERIPIAELF SLDGGSIDAR EDDYEQKAAD PKFTGFHRLE<br>KALFGDNNTK GMDKYADQLY TDVVDLQKRI SELAFPPSKV VGGAAGLIEE<br>VAASKISGEE DRYSHDLWD FQANVEGSQK IVDLLRPQLQ KANPELLAKV<br>DANFKKVDTI LAKYRTKDGf ENYDKLTDAD RNALKGPITA LAEDLAQLRG<br>VLGLD |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | efeO  |
| <b>Protein Names</b>     | Recommended name: Iron uptake system component EfeO   |
| <b>Expression Region</b> | 27-375  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.         |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |