



# Recombinant tRNA 2-thiocytidine biosynthesis protein TtcA (ttcA)

<b>Product Code</b>	CSB-YP372240EJE
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A1AAV5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O1:K1 / APEC
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<pre>MQENQKNTKK EIYNLNKLQK RLRRNVGEAI ADFNMIEEGD RIMVCLSGGK DSYTMLEILR NLQQSAPINF SLVAVNLDQK QPGFPEHVLP EYLETLGVEY KIVEENTYGI VKEKIPEGKT TCSLCSRLRR GILYRTATEL GATKIALGHH RDDILQTLFL NMFYGGKMKG MPPKLMSDDG KHIVIRPLAY CREKDIQRFA DAKAFPIIPC NLCGSQPNLQ RQVIADMLRD WDKRYPGRIE TMFSAMQNVV PSHLCDTNLF DFKGITHGSE VVNGGDLAFD REEIPLQPSG WQPEEDENQL DELRLNVVEV K</pre>
<b>Source</b>	Yeast
<b>Target Names</b>	ttcA
<b>Protein Names</b>	Recommended name: tRNA 2-thiocytidine biosynthesis protein TtcA
<b>Expression Region</b>	1-311
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.