



# Recombinant UDP-N-acetylglucosamine 1-carboxyvinyltransferase (murA)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP372360EJE   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | A1AG91  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O1:K1 / APEC   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MDKFRVQGPT KLQGEVTISG AKNAALPILF AALLAEEPVE IQNVPKLKDV<br>DTSMKLLSQL GAKVERNGSV HIDARDVNVF CAPYDLVKTM RASIWALGPL<br>VARFGQQQVS LPGGCTIGAR PVDLHISGLE QLGATIKLEE GYVKASVDGR<br>LKGAHIVMDK VSVGATVTIM CAATLAEGTT IIENAAREPE IVDTANFLIT<br>LGAKISGQGT DRIVIEGVER LGGGVYRVLP DRIETGTFVLV AAAISRGKII<br>CRNAQPDTLD AVLAKLRDAG ADIEVGEDWI SLDMHGKRPK AVNVRTAPHP<br>AFPTDMQAQF TLLNLVAEGT GFITETVFEN RFMHVPELSR MGAHAEIESN<br>TVICHGVEKL SGAQVMATDL RASASLVLG CIAEGTTVVD RIYHIDRGYE<br>RIEDKLRLALG ANIERVKGE |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | murA  |
| <b>Protein Names</b>     | Recommended name: UDP-N-acetylglucosamine 1-carboxyvinyltransferase<br>EC= 2.5.1.7 Alternative name(s): Enolpyruvate transferase UDP-N-<br>acetylglucosamine enolpyruvyl transferase Short name= EPT  |
| <b>Expression Region</b> | 1-419   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |