



# Recombinant Aspartate carbamoyltransferase (pyrB)

<b>Product Code</b>	CSB-MP372447EJE
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A1AJF0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O1:K1 / APEC
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MANPLYQKHI ISINDLSRDD LNLVLATAAK LKANPQPELL KHKVIASCFF EASTRTRLSF ETSMHR LGAS VVGFS DSANT SLGKKGETLA DTISVISTYV DAIVMRHPQE GAARLATEFS GNV PVLNAGD GSNQHPTQTL LDLFTIQETQ GRLDNLHVAM VGDLKYGRTV HSLTQALAKF DGNRFYFIAP DALAMPQYIL DMLDEKGI AW SLHSSIEEVM AEVDILYMTR VQKERLDPSE YANVKAQFVL RASDLHNAKA NMKVLHPLPR VDEIATDV DK TPHAWYFQQA GNGIFARQAL LALVLNRDLV L
<b>Source</b>	Mammalian cell
<b>Target Names</b>	pyrB
<b>Protein Names</b>	Recommended name: Aspartate carbamoyltransferase EC= 2.1.3.2 Alternative name(s): Aspartate transcarbamylase Short name= ATCase
<b>Expression Region</b>	1-311
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.