



Recombinant Human Clusterin (CLU)

Product Code	CSB-BP005595HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P10909
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	DQTVSDNE LQEMSNQGSK YVNKEIQNAV NGVKQIKTLI EKTNEERKTL LSNLEEAKKK KEDALNETRE SETKLKELPG VCNETMMALW EECKPCLKQT CMKFYARVCR SGSGLVGRQL EEFLNQSSPF YFWMNGDRID SLENDRQQT HMLDVMQDHF SRASSIIDEL FQDRFFTREP QDTYHYLPFS LPHRRPHFFF PKSRIVRSLM PFSPYEPLNF HAMFQPFLEM IHEAQQAMDI HFHSPAFQHP PTEFIREGDD DRTVCREIRH NSTGCLRMKD QCDKCREILS VDCSTNNPSQ AKLRRELDDES LQVAERLTRK YNELLKSYQW KMLNTSSLLE QLNEQFNWVS RLANLTQGED QYYLRVTTVA SHTSDSDVPS GVTEVVVKLF DSDPITVTVP VEVSRKNPKF METVAEKALQ EYRKKHREE
Source	Baculovirus
Target Names	CLU
Protein Names	Recommended name: Clusterin Alternative name(s): Aging-associated gene 4 protein Apolipoprotein J Short name= Apo-J Complement cytolysis inhibitor Short name= CLI Complement-associated protein SP-40,40 Ku70-binding pr
Expression Region	23-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.