



Recombinant Danio rerio Ras-like protein family member 11A-like (zgc:110179)

Product Code	CSB-EP373941DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A1DZY4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MRLLEIQRPT MSSGSSNFLL VPIPEYPVLD CVPNKNVKIV VLGASNVGKT ALIVRFLTKR FIGDYEANTG ALYSRKNLD GEQVSLQVQD TPCVSLQDDA DGLYCQEQIN RSIYWADGYV LVFSITDLNS YRTIQPLYQH VRRIHPSGNI PVIIVGNKSD LLRARQVSDP EGKALADELG GLYFEASARE NHESVQAAFL HLCQEVSRAL GGGNGEKRKG GLHLARPKSP NMQELKRRFR QVLSSKVKSA TAL
Source	E.coli
Target Names	zgc:110179
Protein Names	Recommended name: Ras-like protein family member 11A-like
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.