



# Recombinant Mycobacterium bovis Antigen 85-A (fbpA)

<b>Product Code</b>	CSB-EP374958MOJ-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A1KQD8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mycobacterium bovis (strain BCG / Pasteur 1173P2)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AFSRPGLP VEYLQVSPSPS MGRDIKVQFQ SGGANSPALY LLDGLRAQDD FSGWDINTPA FEWYDQSGLS VVMPVGGQSS FYSDWYQPAC GKAGCQTYKW ETFLTSELPG WLQANRHVKP TGSVVGLSM AASSALTLAI YHPQQFVYAG AMSGLLDPSQ AMGPTLIGLA MGDAGGYKAS DMWGPKEDEPA WQRNDPLLNK GKLIANNTRV WVYCGNGKPS DLGGNNLPAK FLEGFVRTSN IKFQDAYNAG GGHNGVDFDP DSGTHSWEYW GAQLNAMKPD LQRALGATPN TGPAPQGA
<b>Source</b>	E.coli
<b>Target Names</b>	fbpA
<b>Protein Names</b>	Recommended name: Antigen 85-A Alternative name(s): Antigen 85 complex A Short name= 85A Short name= Ag85A Fibronectin-binding protein A Mycolyl transferase 85A EC= 2.3.1.-
<b>Expression Region</b>	43-338
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.