



Recombinant Danio rerio Serine/threonine-protein phosphatase 4 regulatory subunit 2-B (ppp4r2b)

Product Code	CSB-MP375502DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A1L1M4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MEIDTLLEAL KDFDKKVKKE ASTVLDQFLC HIAKTGETLV SWSQFKTYFV FKLEKVMEDF HASAPEQRGP ANPNVESIPF EEMKERILKI VNEYNGIPFT IQRLCELLTE PKRNYTGTEK FLRGVEKNVM VVSCIYPTSE KNGCSGVNRM NGVMLSGNTS AFTERKVNGP GTPRPLNRQK HSMSSSLATN GLPDSTENKD TNT PQGHSKT TSDASASGIG GPLGSSVKNK HLDDEEDEGM EEEEHEVKRL KFSQDGNED EDEDDVDEQV VETPRSTLSK SPSPDKSEDPE SPSSMLQEGS EEEEESKHE GSAENSEDQE PSSTQPECPV ADSRETFTER EVPCGSDSS EEPSDMDQSE QRAPAGVASS PELREGEEPS EVSSSSTDDN ASGPTSSSSS IDSPTEGDTG SGALDLESAE EPMEQD
Source	Mammalian cell
Target Names	ppp4r2b
Protein Names	Recommended name: Serine/threonine-protein phosphatase 4 regulatory subunit 2-B
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.