



Recombinant Arabidopsis thaliana Alcohol dehydrogenase-like 3 (At1g32780)

Product Code	CSB-YP375946DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A1L4Y2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAETQGKVVIT CKAADVWVWVWPK VPLVIQVQICV DPPQKMEVVRV KILYSSICHT DLGCWNGTNE AERAFPRILG HEAVGIVESV GEGVKDVKEG DYVIPTFNGE CGECKVCKRE ESNLCERYHV DPMKRVMVND GGTRFSTTIN KGGSSQSQP IYHFLNTSTF TEYTVLDSAC VVKIDPNSPL KQMSLLSCGV STGVGAAWNI ANVKEGKSTA VFGLGSVGLA VAEGARARGA SRIIGVDANA SKFEKGKLMG VTDFINPKDL TKPVHQMIRE ITGGVDYSF ECTGNVDVLR EAFSTHVGW GSTVLVGIYP TPRTLPLHPM ELFDGRRITG SVFGGFKPKS QLPNFAQQCM KGVVKLEPFI TNELPFKIN DAFQLLRDVGK SLRCILQISK LLKR
Source	Yeast
Target Names	At1g32780
Protein Names	Recommended name: Alcohol dehydrogenase-like 3 EC= 1.1.1.1
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.