



Recombinant Chicken Proto-oncogene c-Fos (FOS)

Product Code	CSB-EP008790CH-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P11939
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	<p>MMYQGFAGEY EAPSSRCSSA SPAGDSLTY Y PSPADSFSSM GSPVNSQDFC TDLAVSSANF VPTVTAISTS PDLQWLVQPT LISSVAPSQN RGHPYGVPA P APPAAYS RPA VLKAPGGRGQ SIGRRGKVEQ LSPEEEEEKRR IRRERNKMAA AKCRNRREL TDTLQAETDQ LEEEKSA LQA EIANLLKEKE KLEFILAAHR PACKMPEELR FSEELAAATA LDLGAPSPAA AEEAFALPLM TEAPPAVPPK EPSGSGLELK AEPFDELLFS AGPREASRSV PDMDLPGASS FYASDWEPLG AGSGGELEPL CTPVVTCTPC PSTYTSTFVF TYPEADAFPS CAAHRKGS S SNEPSSDSL S SPTLLAL</p>
Source	E.coli
Target Names	FOS
Protein Names	Recommended name: Proto-oncogene c-Fos Alternative name(s): Cellular oncogene fos
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differentiation, and transformation. In some cases, expression of the FOS gene has also been associated with apoptotic cell death.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.