



Recombinant Human Teratocarcinoma-derived growth factor 1 (TDGF1)

Product Code	CSB-EP023343HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P13385
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	LGHQEFARPS RGYLAFRDDS IWPQEEP AIR PRSSQRVPPM GIQHSKELNR TCCLNGGTCM LGSFCACPPS FYGRNCEHDV RKENCGSVPH DTWLPPKKCSL CKCWHGQLRC FPQAFLPGCD
Source	E.coli
Target Names	TDGF1
Protein Names	Recommended name: Teratocarcinoma-derived growth factor 1 Alternative name(s): Cripto-1 growth factor Short name= CRGF Epidermal growth factor-like cripto protein CR1
Expression Region	31-150
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.