



Recombinant Danio rerio Uncharacterized protein FLJ46082 homolog (si:ch211-215c18.4)

Product Code	CSB-BP381198DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2BFC9
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MESPRIGVVR DSMLHRPLLI KPTLGKSKSR GLSYPGPEFV FGTATTVQDG GVQEAISNWH THTMSTTHRA AERDFIALNR EGVKSGLVTA KELHQYRATH DIRRKPLTRE GFRRSAPPRL SPDDSFVSN RPSTPISELI EYKYAQRWLE EQQVRDRVLQ AHQHKKAQLG RIQDTRTLL RKSQPLPEAP STWKLPRFQK VGAALDTFRD PEARKKAMNA HCSESVSRRG ILGQGTYTVD
Source	Baculovirus
Target Names	cfap77
Protein Names	Recommended name: Uncharacterized protein FLJ46082 homolog
Expression Region	1-240
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.