



Recombinant Danio rerio Serologically defined colon cancer antigen 3 homolog (sdccag3)

Product Code	CSB-EP381498DIL-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2CE83
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSRHKSKKLI IEDDEPKQDE LNPFSFKEFI RNKNQPSCST EATEVYNGAC LVEQDYNTSI DFATKGPFFT DPSALCQPPE SEPDETWIES HQRSAIEGSP DFELCGASDI SAYSEQSSLC SDEREAAEWE LGQSDFLPKK HLSRISTGSY EGDEETSVID ISFHHRGNA ENGTKNLQQL REENSLLRKQ VKELVRMSET DSRRIKQLTD EMHNKKLQEE KEANDLEAMV QSVEQNLQLM TKRAVKAENN ISKCLKQEMLK LQGQLEDYKS ENERLRLGET AALATMRHNA QVASEYLNKA AQDAEISINS C
Source	E.coli
Target Names	sdccag3
Protein Names	Recommended name: Serologically defined colon cancer antigen 3 homolog
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.