



# Recombinant Human Keratin, type I cytoskeletal 13 (KRT13)

<b>Product Code</b>	CSB-MP012511HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P13646
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSLRLQSSSA SYGGGFGGGS CQLGGGRGVS TCSTRFVSGG SAGGYGGGVS CGFGGGAGSG FGGGYGGGLG GGYGGGLGGG FGGGFAGGFV DFGACDGGLL TGNEKITMQN LNDRLASYLE KVRALEEANA DLEVKIRDWH LKQSPASPER DYSPYYKTIE ELRDKILTAT IENNRVILEI DNARLAADD F RLKYENELAL RQSVEADING LRRVLDELTL SKTDLEMQIE SLNEELAYMK KNHEEEMKEF SNQVVGQVNV EMDATPGIDL TRVLAEMREQ YEAMAERNRR DAEWFHTKS AELNKEVSTN TAMIQTSKTE ITELRRTLQG LEIELQSLS MKAGLENTVA ETECRYALQL QQIQGLISSI EAQLSELRSE MECQNQEYKM LLDIKTRLEQ EIATYRSLE GQDAKMIGFP SSAGSVSPRS TSVTTTSSAS VTTTASNASGR RTSDVRRP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	KRT13
<b>Protein Names</b>	Recommended name: Keratin, type I cytoskeletal 13 Alternative name(s): Cytokeratin-13 Short name= CK-13 Keratin-13 Short name= K13
<b>Expression Region</b>	1-458
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This protein is a member of the keratin gene family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. Most of the type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. This type I cytokeratin is paired with keratin 4 and expressed in the suprabasal layers of non-cornified stratified epithelia. Mutations in this gene and keratin 4 have been associated with the autosomal dominant disorder White Sponge Nevus. The type I cytokeratins are clustered in a region of chromosome 17q21.2. Alternative splicing of this gene results in multiple transcript variants; however, not all variants have been described.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.