



Recombinant Danio rerio UPF0472 protein C16orf72 homolog (si:dkey-20i6.3, zgc:77849)

Product Code	CSB-BP382950DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2BE76
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MEEKKEEGEA EIQEHGPEHW FSKWERQCLA EAEQEEPAAE ETDQSQQKLW HLFQNSATAV AQLYKDRVCQ QQQGLSLWVP FQNAATAVTN LYKESVDAHQ RSYELGIQIG HQRRNKDVLA WVKKRRRTIR REDLISFLCG KAPPPRSSRA PPRLAMVSPS RSTPSETSSS VETDLQPFRE AIALHGLSGA MASISMRS GA PGSPTHLSAS SAPSRRRNGL HDVDLNTFIA EEMALHLDNG TRKRSSAQCN DVITDSPTHK RNRMI
Source	Baculovirus
Target Names	si:dkey-20i6.3
Protein Names	Recommended name: UPF0472 protein C16orf72 homolog
Expression Region	1-265
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.