



Recombinant *Oryza sativa* subsp. *indica* Dehydration-responsive element-binding protein 2A (DREB2A)

Product Code	CSB-YP383680OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2WL19
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MERGEGRGD CSVQVRKKRT RRKSDGPDSI AETIKWWKEQ NQKLQEENSS RKAPAKGSKK GCMAGKGGPE NSNCAYRGVR QRTWGKWVAE IREPNRGRRL WLGSFPTALE AAHAYDEAAR AMYGP TARVN FADNSTDANS GCTSAPSLMM SNGPATIPSD EKDELESPPF IVANGPAVLY QPDKKDVLER VVPEVQDVKT EGSNGLKRVC QERKTMEVCE SEGIVLHKEV NISYDYFNVH EVVEMIIVEL SADQKTEVHE EYQEGDDGFS LFSY
Source	Yeast
Target Names	DREB2A
Protein Names	Recommended name: Dehydration-responsive element-binding protein 2A Short name= Protein DREB2A
Expression Region	1-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.