



Recombinant *Oryza sativa* subsp. *indica* COBRA-like protein 5 (BC1)

| | |
|--------------------------|--|
| Product Code | CSB-BP383706OFF |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | A2XHZ9 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Oryza sativa</i> subsp. <i>indica</i> (Rice) |
| Purity | >85% (SDS-PAGE) |
| Sequence | YDPLDPKG NITIKWDVIS WTPDGYVAMV TMSNYQMYRQ ILAPGWTVGW SWAKKEVIWS IVGAQATEQG DCSKFKGGIP HCKRTPAIV DLLPGVPYNQ QIANCCKAGV VSAYGQDPAG SVSAFQVSVG LAGTTNKTVK LPTNFTLAGP GPGYT�CGPAT IVPSTVYLTP DRRRRTQALM TWTVTCTYSQ QLASRYPTCC VSFSSFYNST IVPCARCACG CGHDGYRGNG GGGKNARAGD GRSRNSGGG GGHSGGTECI MGDSKRALSA GVNTPRKDGA PLLQCTSHMC PIRVHWHVKL NYKDYWRAKI AITNFNYRMN YQWTLVAQH PNLNNVTEVF SFQYKPLLPY GNINDTGMFY GLKFYNDLLM EAGPFGNVQS EVLMRKDYNT FTFSQGWAFP RKIYFNGDEC KMPPPSYPY LPN |
| Source | Baculovirus |
| Target Names | BC1 |
| Protein Names | Recommended name: COBRA-like protein 5 Alternative name(s): Protein BRITTLE CULM1 |
| Expression Region | 23-443 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |