



Recombinant *Oryza sativa* subsp. *indica* COBRA-like protein 5 (BC1)

Product Code	CSB-YP383706OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2XHZ9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	YDPLDPKG NITIKWDVIS WTPDGYVAMV TMSNYQMYRQ ILAPGWTVGW SWAKKEVIWS IVGAQATEQG DCSKFKGGIP HCKRTPAIV DLLPGVPYNQ QIANCCKAGV VSAYGQDPAG SVSAFQVSVG LAGTTNKTVK LPTNFTLAGP GPGYT�CGPAT IVPSTVYLTP DRRRRTQALM TWTVTCTYSQ QLASRYPTCC VSFSSFYNST IVPCARCACG CGHDGYRGNG GGGKNARAGD GRSRNSGGG GGHSGGTECI MGDSKRALSA GVNTPRKDGA PLLQCTSHMC PIRVHWHVKL NYKDYWRAKI AITNFNYRMN YQWTLVAQH PNLNNVTEVF SFQYKPLLPY GNINDTGMFY GLKFYNDLLM EAGPFGNVQS EVLMRKDYNT FTFSQGWAFP RKIYFNGDEC KMPPPSYPY LPN
Source	Yeast
Target Names	BC1
Protein Names	Recommended name: COBRA-like protein 5 Alternative name(s): Protein BRITTLE CULM1
Expression Region	23-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.