



# Recombinant *Oryza sativa* subsp. *indica* Actin-3 (ACT3)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP383713OFF                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                  |
| <b>Uniprot No.</b>       | A2XNS1                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Immunogen Species</b> | <i>Oryza sativa</i> subsp. <i>indica</i> (Rice)                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Sequence</b>          | MADGEDIQPL VCDNGTGMVK AGFAGDDAPR AVFPSIVGRP<br>RHTGVMVGMG QKDAYVGDEA QSKRGILTLK YPIEHGIVNN WDDMEKIWHH<br>TFYNELRVAP EEHPILLTEA PLNPKANREK MTQIMFETFN APAMYVAIQA<br>VLSLYASGRT TGIVLDSGDG VTHTVPIYEG YALPHAILRL DLAGRDLTDC<br>LMKILTERGY SFTTTAEREI VRDIKEKLAY IALDYEQELE TAKSSSSVEK<br>SYELPDGQVI TIGAERFRCP EVLFQPSLIG MEAPGIHETT YNSIMKCDVD<br>IRKDLYGNIV LSGGTTMFPG IADRMSKEIT ALAPSSMKIK VVAPPERKYS<br>VWIGGSILAS LSTFQQMWIS KGEYDESGPS IVHRKCF |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Target Names</b>      | ACT3                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Protein Names</b>     | Recommended name: Actin-3                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Expression Region</b> | 1-377                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                 |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                         |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                  |