



Recombinant *Oryza sativa* subsp. *indica* Probable protein arginine N-methyltransferase 6.2 (PRMT6.2)

Product Code	CSB-BP383755OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2Z8S0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MFAGGADGGN GHLPRRRAR RGGGGGGGMG SPPLGPPPPP CTDYDMAYFK AYSHIGVHEE MLKDHVRTNT YRNAIMHHQD LISGKVLDV GCGTGVSIF CAFAGAARVY AVDASDIALQ AMEIVRENEL SDKVIVLHGR IEDVEIEEKV DVIISEWMGY MLLYESMLGS VIFARDKWLK PGGLILPSHA SLYLAPITNS HRYQDSVYFW QDVYGIKMSS MMPLAKQCAF MEPSVETISG ENVLTWPSVV AQVDCYTIQA PELETITATF NYTSMLQAPL HGFAFWFDVE FNGPVRQRSK KQANQCLDGN TQDASPSNKK KKADAPIVLS TAPEDAPTHW QQTLLYLFEP IELKKDQNE GSVTISQSQQ HARFLNICKL YFTRDQWYVK ESVMK
Source	Baculovirus
Target Names	PRMT6.2
Protein Names	Recommended name: Probable protein arginine N-methyltransferase 6.2 EC=2.1.1.-
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.