



Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 9 (Os01g0645000, LOC_Os01g45730)

Product Code	CSB-MP383765OFG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2ZVY5
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MQEALVSPVR SALELKPFNF GDQRLASSPR YLPSGDDALY RCSSPFSPSF GFSSPSPLAT SVSLSPSSSA SLVDDGDDGG AAADATGQRL QLARLALQYQ EVADRYELCL SHLAEAAEEA AALRLENAEL RVTNSDLALR LALLSGKHTA AVAVADEIRR LRLGEQKVA AATKERTPEKL AVLPKSISVR STSYLKLNQQ SQAATATSAA PNRKPRTSSN PTNPPNSQRA YDGGKKGDEQ KAAPADSGAE LEVYNQGMFK TELCNKWEET GDCPYGDQCQ FAHGVTELRP VIRHPRYKTA VCRMVLAGDV CPYGHRCHFR HSLTPAERLL LRS
Source	Mammalian cell
Target Names	Os01g0645000
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 9 Short name= OsC3H9
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.