



Recombinant Scheffersomyces stipitis Protein transport protein SEC13 (SEC13)

Product Code	CSB-MP384103SXM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A3LNW3
Product Type	Recombinant Protein
Immunogen Species	Scheffersomyces stipitis (strain ATCC 58785 / CBS 6054 / NBRC 10063 / NRRL Y-11545) (Yeast) (Pichia stipitis)
Purity	>85% (SDS-PAGE)
Sequence	MVTIGNAHND LIHDAVLDDYY GKRLATCSSD KSINIFDIDG TESYKLVSTL TGHGDPVWQV SWAHPKFGSI LASCSEFDGKA LIWKEQPETQ QWSIIAEHSV HQASVNSVSW APHELGA VLL CASSEGDGKSV VDFNDDGTTS HVVFDHAHAIG ANSASWAPLS STPSPNQKDA ASLKQQRFFV TCGSDNLAKI WKYDAANNTY VEEARLEGHT DWVRDVAWSP SMLVRTYIAT ASQDRTVLIW TQDKAGKWQK QLLTEDKFPD VCWRCSWSLS GNILAVSGGD NKVSLWKENL QGKWESAGEV VQ
Source	Mammalian cell
Target Names	SEC13
Protein Names	Recommended name: Protein transport protein SEC13
Expression Region	1-302
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.