



# Recombinant Callithrix jacchus Glutathione reductase, mitochondrial (GSR)

<b>Product Code</b>	CSB-BP384847CYL
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A2TIL1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Callithrix jacchus (White-tufted-ear marmoset)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MACRQEP QPQGPPPSAG AVVSYDYLVI GGGSGGLASA RRAAELGARA AVVESHKLGG TCVNVGCVPK KVMWNTAVHS EFLHDHGDYG FSSCEGKFNW RVIKEKRDTY VSRLNTIYQN NLTKAHIEII HGHAVFTSDT KPTIEVSGRK YTAPHILIAT GGMPSSPHES QIPGASLGIT SDGFFELEEL PSRSVIVGAG YIAVEIAGIL SALGSKTSLM IRHDKVLRSF DSMISTNCTE ELENAGVEVL KFSQVKEVKK TSSGLEVSLV TAVPGRLPVM TTISDVDCLL WAIGRDPNSK GLSLNKLGIK TDDKGHIIVD EFQNTNVKGI YAVGDVCGKA LLTPVAIAAG RKLAHRLFEN KEDSKLDYNN IPTVVFSHPP IGTVGLTEDE AIHKYGKENV KIYSTSFTPM YHAVTKRRTK CVMKMVCAYE EEKVVGIHMQ GLGCDEMLQG FAVAVKMGAT KADFDNTVAI HPTSSEELVP LR
<b>Source</b>	Baculovirus
<b>Target Names</b>	GSR
<b>Protein Names</b>	Recommended name: Glutathione reductase, mitochondrial Short name= GR Short name= GRase EC= 1.8.1.7
<b>Expression Region</b>	44-522
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.