



Recombinant *Oryza sativa* subsp. *indica* Uncharacterized protein Osl_009113 (Osl_009113)

Product Code	CSB-EP384871OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2XAM0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MAQNKTIAVA LLLATLVAVM GKEPETLEET LragckEECS EQKKKAPIDE KQCEDFCFIK TKSIFEAHKG VKDLKADRFI DFCNNECNAV YKEDPATSKK CAESCEADAK EAEVFLDKVV AYMQTMKQA
Source	<i>E.coli</i>
Target Names	Osl_009113
Protein Names	Recommended name: Uncharacterized protein Osl_009113 Alternative name(s): Unknown protein AN04 from 2D-PAGE of anther
Expression Region	1-129
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.