



Recombinant *Oryza sativa* subsp. *indica* Auxin-responsive protein IAA10 (IAA10)

Product Code	CSB-MP384872OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2XB18
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MRGGAAGPTA GEPGTEAEA EEVEESSAGD DEEELGLSL GSKKQQQQQH APCRILTARD LQPAAALSPD SSVSSSSPAA AAAAGGKRAE GPTATTSPGT VASGHPHSSF GVVGWPPIRQ FRMNSLNFQA KENTSETDTK KTATNESDVQ KDKEEGEKKG RVAGWVKVNM DGEVIGRKVD LNAHRSYKTL ALALELMFTK PSIGLCASHN TNSLKLDDNS AEYQLTYEDR DGDWMLVGDV PWEMFVSSVK RLRIMRTSDA NGLGQRYQGI HRTIASTRGR S
Source	Mammalian cell
Target Names	IAA10
Protein Names	Recommended name: Auxin-responsive protein IAA10 Alternative name(s): Indoleacetic acid-induced protein 10
Expression Region	1-281
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.