



# Recombinant *Oryza sativa* subsp. *indica* Mitogen-activated protein kinase 5 (MPK5)

<b>Product Code</b>	CSB-MP384876OFF
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A2XFC8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDGAPVAEFR PTMTHGGRYL LYDIFGNKFE VTNKYQPPIM PIGRGAYGIV CSVMNFETRE MVAIKKIANA FNNDMDAKRT LREIKLLRHL DHENIIGIRD VIPPPIPQAF NDVYIATELM DTDLHHIIRS NQELSEEHCQ YFLYQILRGL KYIHSANVIH RDLKPSNLLL NANCCLKICD FGLARPSSES DMMTEYVVTR WYRAPELLLN STDYSAAIDV WSVGCFMEL INRQPLFPGR DHMHQMRLIT EVIGTPTDDE LGFIRNEDAR KYMRHLPQYP RRTFASMFPR VQPAALDLIE RMLTFNPLQR ITVEEALDHP YLERLHDIAD EPICLEPFSF DFEQKALNED QMKQLIFNEA IEMNPNIRY
<b>Source</b>	Mammalian cell
<b>Target Names</b>	MPK5
<b>Protein Names</b>	Recommended name: Mitogen-activated protein kinase 5 Short name= MAP kinase 5 EC= 2.7.11.24 Alternative name(s): Benzothiadiazole-induced MAP kinase 1 MAP kinase 2 Multiple stress-responsive MAP kinase 2 OsBIMK1 OsMAP
<b>Expression Region</b>	1-369
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.