



Recombinant *Oryza sativa* subsp. *indica* WUSCHEL-related homeobox 6 (WOX6)

Product Code	CSB-MP384878OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2XG77
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEGSSNSPDR QSSGGSPPEE RGGGGSGGGG GRSAAGEPVR SRWTPKPEQI LILESIFNSG MVNPPKDETV RIRKLLERFG AVGDANVFYW FQNRRSRSRR RQRQMMAAAA AAAAASSSS PSANTSPAAA SAATVQVGLP PVAVVHTMAM GGSACQYEQQ ASSSSSSGST GGSSLGLFAH GAGASGAGGY LQASCGASAS ASSALAPGLM GDVVDSSGGSD DLFAISRQMG FVGSPRCSPA SSPATPSSAA TAAQQQFYSC QLPAATITVF INGVPMEMPR GPIDLRAMFG QDVMLVHSTG ALLPVNDYGI LMQSLQIGES YFLVARPT
Source	Mammalian cell
Target Names	WOX6
Protein Names	Recommended name: WUSCHEL-related homeobox 6 Alternative name(s): OsWOX6
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.