



Recombinant *Oryza sativa* subsp. *indica* Transcription factor PCF6 (PCF6)

Product Code	CSB-YP384884OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2XMN1
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEAAVGDGEG GGGGGGRGKR GRGGGGGEMV EAVWGQTGST ASRIYRVRAT G GKDRHSKVY TAKGIRD RR V RLSVATAIQF YDLQDRLGFD QPSKAIEWLI NAASPAIDTL PSLDPAAF AA IPHAAAADAA PTRRRSQQQQ QQLSNKSGCS STSETSKGSD KEVTVASAPA QAASFTELLI AGVAASSAGG GAINGADCV GIAHPGKGG A EGASTYG FSA ASSFGDAPPI GMVPAPPFNF SAPGADMAAH YSLAQDQLAA PPPPAGGDYN LNFSMSSGFL GANRGTLQSN SPSNMSGHHH HHHQQQLQRL DGSTISFLG HAAAAHPAA SEGQITSTAA LQLWDGFRHS GMKEKSKN
Source	Yeast
Target Names	PCF6
Protein Names	Recommended name: Transcription factor PCF6
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.