



Recombinant *Oryza sativa* subsp. *indica* S-adenosylmethionine synthase 1 (SAM1)

Product Code	CSB-BP384893OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2Y053
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAALDTFLFT SESVNEGHPD KLCDQVSDAV LDACLAEDPD SKVACETCTK TNMVMVFGEI TTKANVDYK IVRETCRNIG FVSADVGLDA DHCKVLVNI QQSPDIAQGV HGHFTKRPEE IGAGDQGHMF GYATDETPEL MPLSHVLATK LGARLTEVRK NGTCAWLRPD GKTQVTVEYR NESGARVPVR VHTVLISTQH DETVTNDIEIA ADLKEHVIK VPEQYLDEK TIFHLNPSGR FVIGGPHGDA GLTGRKIID TYGGWGAHGG GAFSGKDPTK VDRSGAYVAR QAAKSIVASG LARRCIVQVS YAIGVPEPLS VFDVTYGTGR IPDKEILKIV KENDFRPGM IINLDLKKG GNGRYLKTA YGHFGRDDPD FTWEVVKPLK WEKPSA
Source	Baculovirus
Target Names	SAM1
Protein Names	Recommended name: S-adenosylmethionine synthase 1 Short name= AdoMet synthase 1 EC= 2.5.1.6 Alternative name(s): Methionine adenosyltransferase 1 Short name= MAT 1
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.