



Recombinant *Oryza sativa* subsp. *indica* Probable acetyltransferase NSI (NSI)

Product Code	CSB-EP384897OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2Y5T7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MASAASASAS AVVTPSSFRC VPTASCGLGA RGKAPAPRRL LHDHAQGKKR AAATWSLKAG LWDSLRSGLF KSNSTETVE PPSAPIEEEE PLPEELVLE RTLADGSTEQ IIFSSAGDVN VYDLQALCDK VGWPRRPLTK IAASLRNSYL VATLHSV TMP SKAEGEERKQ LIGMARATSD HAFNATIWDV LVDPSYQQGG LGKALMEKVI RTLLQRDISN ITLFADNKVV DFYKNLGF EA DPQGIKGMFW YPRF
Source	<i>E. coli</i>
Target Names	SNAT1
Protein Names	Recommended name: Probable acetyltransferase NSI EC= 2.3.1.- Alternative name(s): Nuclear shuttle protein-interacting protein homolog
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.