



Recombinant *Oryza sativa* subsp. *indica* Catalase isozyme B (CATB)

Product Code	CSB-MP384906OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2YH64
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDPYKHRPSS GSNSTFWTTN SGAPVWNNNS ALTVGERGPI LLEDYHLIEK LAQFDRERIP ERVVHARGAS AKGFFEVTHTD ISHLTCADFL RAPGVQTPVI VRFSTVVHER GSPETLRDPR GFAVKFYTRE GNFDLVGNM PVFFIRDGMK FPDMVHAFKP SPKTNMQENW RIVDFFSHHP ESLHMFSLF DDVGIPLNYR HMEGFGVNTY TLINKDGKPH LVKFHWKPTC GVKCLLDDEA VTVGGTCHSH ATKDLTDSIA AGNYPEWKLY IQTIDPDHED RFDLPLDVT KTWPEDIPL QPVGRMVLNK NIDNFFAENE QLAFCPAIIV PGIHYSDDKL LQTRIFSYAD TQRHRLGPNY LMLPVNAPKC AYHNNHHDGS MNFMHRDEEV NYFPSRFDAA RHAEKVPIPP RVLTCREKC VIDKENNFKQ AGERYRSFDP ARQDRFLQRW VDALS DPRIT HELRGIWISY WSQCASLGQ KLASRLNLKP NM
Source	Mammalian cell
Target Names	CATB
Protein Names	Recommended name: Catalase isozyme B Short name= CAT-B EC= 1.11.1.6
Expression Region	1-492
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.