



Recombinant *Oryza sativa* subsp. *indica* Serine/threonine-protein kinase SAPK2 (SAPK2)

Product Code	CSB-EP384911OFF-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2YNT8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MERYEVIKDI GSGNFGVAKL VRDVRTKELF AVKFIERGQK IDENVQREIM NHRSLRHPNI VRFKEVVLTP THLAIVMEYA AGGELFERIC SAGRFSEDEA RFFFQQLISG VSYCHSMQIC HRDLKLENTL LDGSIAPRLK ICDFGYSKSS VLHSQPKSTV GTPAYIAPEV LARKEYDGKV ADVWSCGVTL YVMLVGAYPF EDPDEPRNFR KTITRILSVQ YMVPDYVRVS MECRHLLSRI FVANPEQRIT IPEIKNHPWF LKNLPIEMTD EYQMSVQMND INTPSQGLEE IMAIQEARL PGDGSKFSGQ IPGLGSMELD DIDTDDIDVE DSGDFVCAL
Source	E.coli
Target Names	SAPK2
Protein Names	Recommended name: Serine/threonine-protein kinase SAPK2 EC= 2.7.11.1 Alternative name(s): Osmotic stress/abscisic acid-activated protein kinase 2
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.