



Recombinant *Oryza sativa* subsp. *indica* Zinc finger A20 and AN1 domain-containing stress-associated protein 1 (SAP1)

Product Code	CSB-MP384924OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2Z2J6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAQRDKKDQE PTELRAPEIT LCANSCGFPG NPATQNLCQN CFLAATASTS SPSSLSSPVL DKQPPRPAAP LVEPQAPLPP PVEEMASALA TAPAPVAKTS AVNRCSRKRK RVGLTGFRRCR CGHLFCGEHR YSDRHGCSYD YKSAARDAIA RDNPVVRAAK IVRF
Source	Mammalian cell
Target Names	SAP1
Protein Names	Recommended name: Zinc finger A20 and AN1 domain-containing stress-associated protein 1 Short name= OsSAP1 Alternative name(s): Multiple stress-responsive zinc finger protein ISAP1 Short name= OsiSAP1
Expression Region	1-164
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.