



Recombinant *Oryza sativa* subsp. japonica

Putative protein phosphatase 2C 76

(Os11g0586001, LOC_Os11g37540)

Product Code	CSB-BP384950OFG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A3CCP9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	ESATCL AVYREGGAPA VFQSAHCPRW TLLAPSAGSG GEGDGDRRSS SSSPPPPPHP RGCHVAVDRG RRRSQEDRAV CALGIRIPFI GIYHKIKEVD VGVVAVFDGH NGAEASEMAS KLLLEYFLLH VYFLLDGIYS IMFRKSTGKL TYKEVTILNN VINLYKEDQS SHSKGSCWAL PAILDRSFHM EVLKESLLRA VHDVDLTFSK EALRNNEFESG STAAVILIVD GQIIAANVGD SKAFLCSESH DSSIDKKTSV VSGKRRRKRN SNNRDDFALA NYDGPFFYNVK ELTKDHHPDR EDERSRVEAA GGYVLEWAGV HRVNGELALS RAIGDVPYKR YGVIPTPELT EWQSLSANDT FLIASSDGVF EKMTMQDVCD LMLRVKLGVN QELGSFAVTQ RNLADYVVDL ALEKGTDDNV AAVIVPLGSH YSSKVTLEDW YMLEENSKTS ISPLQTIPYQ QKSGRFNVQ
Source	Baculovirus
Target Names	Os11g0586001
Protein Names	Recommended name: Putative protein phosphatase 2C 76 Short name= OsPP2C76 EC= 3.1.3.16
Expression Region	35-499
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.