



Recombinant Danio rerio Probable methyltransferase C20orf7 homolog, mitochondrial (zgc:162919)

Product Code	CSB-BP385277DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A3KP37
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	VFDRS MKRRQKDWAS SLLDSSKYDY LREEVGSRVA DRVYDVARTF PLALDVGCGR SHIAEHSKE VVERLFLTDI SSSSLRNRKT SDIPAQCVMMA DEEFLPFKEN TFDLVLSSLS MHWINDLPGA LRQIHQVLKP DGVFIGAMVG GETLYELRCS LQLAELEREG GFAPHISPYT AVTDLGNLLG QAGFNMLTVD IDEVQVNYPG MLEVMRDLQG MGESNCAWNR KLLLQRDTML AAAAIYKEMY GNEDGSVPAT FQILYMIGWK PHDSQAKPAK RGSANVSFAD LSKIGKLQSD Q
Source	Baculovirus
Target Names	ndufaf5
Protein Names	Recommended name: Probable methyltransferase C20orf7 homolog, mitochondrial EC= 2.1.1.-
Expression Region	26-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.