



Recombinant *Oryza sativa* subsp. *indica* Metallothionein-like protein 3A (MT3A)

Product Code	CSB-YP386051OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2WLS0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSDKCGNDCD ADKSQCVKKG TSYGVVIVEA EKSHFEEVAA GEENGGCKCG TSCSCTDCKC GK
Source	Yeast
Target Names	MT3A
Protein Names	Recommended name: Metallothionein-like protein 3A Alternative name(s): Class I metallothionein-like protein 3A OsMT-I-3a Short name= OsMT3 Short name= OsMT3a
Expression Region	1-62
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.