



Recombinant *Oryza sativa* subsp. *indica* Squamosa promoter-binding-like protein 3 (SPL3)

Product Code	CSB-MP386063OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2X0Q6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	<p> MGSFGMDWNQ KSSVLWDWEN MPPIGNSANE NPKNVMLAES KLAGVGVDIG HESGHSSGGT FSSSSEIGYG SSKSSISASI DSPSKVGNTI ELNFASAEH DKNMDKGKSK VDDTGTSRSP VVAANRVEPL IGLKLGKRTY FEDVCGGQNV KSSPSGVSVV TPSPGLAKKV KVAQQNTQNP HCQVEGCNVD LSSAKPYHRK HRVCEPHSKT LKVIVAGLER RFCQQCSRFB GLAEFDQKKR SCRRRLHDHN ARRRKPQPEA ISLSSRLST LLYGDARQQA SFLFGQAPYG QMGSCASSWD NPVPGGFKFT ATKAPWSRPT IAAGVDGTHV SNQQASGNVL PHGAHHSFDG LMAFKETNAK VLNQGMEASA VASGSARGPD FEHALSLLSI DSVGAANLQP GSQIHPGVTA IAGTSNPVMM PSPAIWQGGL SLDQQAQFQA FDRLGNDDDE DHLQLPKPSY DNSHYDQMN </p>
Source	Mammalian cell
Target Names	SPL3
Protein Names	Recommended name: Squamosa promoter-binding-like protein 3
Expression Region	1-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.