



Recombinant *Oryza sativa* subsp. *indica* Homeobox-leucine zipper protein HOX13 (HOX13)

Product Code	CSB-MP386073OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2XDD6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MKRPTSSSRK SKKQGEDLAF SEEGSLPAVT MEQKDEAEME EVDDEEEEEV DEDMAGGHAA QSPSPSCGLG EKKRRLALEQ VRALERSFDT DNKLDPDRKA RIARDLGLQP RQVAVWFQNR RARWKTKQLE RDFAALRAQH NDALRADCD LRRDKDALAA EIRELREKLP TKPADTAASV KVEAGNDAAA GAAAATVCKD GSSDDSDSSV VFNDEASPYS GAAFIGFGPS FLVDDASAAT VGCSSSLPAL ESKWHGPYSD DSCKGGVYGF TEEWLAACSG EMAGNDAAGF FSDEHASNLN FGWCASGNEG WE
Source	Mammalian cell
Target Names	HOX13
Protein Names	Recommended name: Homeobox-leucine zipper protein HOX13 Alternative name(s): HD-ZIP protein HOX13 Homeodomain transcription factor HOX13 OsHox13
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.