



# Recombinant *Oryza sativa* subsp. *indica* Probable aldo-keto reductase 2 (Osl\_15387, OSIGBa0115A19.4)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-BP386084OFF  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | A2XRZO   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | <i>Oryza sativa</i> subsp. <i>indica</i> (Rice)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MAAAPATAA VRRMKLGSQG LEVSAQGLGC MGMSAFYGGP KPEPDMVALI<br>HHAVAAGVTL LDTSDIYGPH TNELLGKAL QGGVRDKVEL ATKFGIAFED<br>GKRGVRGDPA YVRAACEGSL RRLGVDSIDL YYQHRVDKKV PIEVTIGELK<br>KLVEEGKIKY IGLSEASAST IRRHAVHPI TAVQLEWSLW SRDVEEDIIP<br>TCRELIGIV AYSPLGRGFF SAGAKLVESL SDQDFRKHIP RFQQENLEKN<br>AEIFERNAM AARKGCTPSQ LALAWVHHQG SDVCPPIGTT KIENLNQNIQ<br>ALSVKLTPEE MAELESYAST DDVRGDRYPQ AMANTTWQNS ETPPLSSWKA<br>Q |
| <b>Source</b>            | Baculovirus  |
| <b>Target Names</b>      | Osl_15387  |
| <b>Protein Names</b>     | Recommended name: Probable aldo-keto reductase 2 EC= 1.1.1.-   |
| <b>Expression Region</b> | 1-351  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.        |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |