



Recombinant *Oryza sativa* subsp. *indica* Anamorsin homolog 1 (H0403D02.15, OsI_17896)

Product Code	CSB-BP386089OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2XYW4
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAATAAAALA VTDELALPLR AVGDAAAAG VSREEVVIT QCASLGGKLP FADASVGSVL AVIKKVENLG NQFITEISRV LKAGGMVLVQ SSPSDQDPNN SIERKLLLGG FVDVQASAAS SQDNEHSVNI KAKKASWSMG SSFPLKKATK GLPKIQIDDD SELIDEDSLL TEDDLKKPEL PVVGDCEVGA TRKACKNCTC GRAEAEKVE KLNLTSEQIN NPQSACGNCG LGDAFRCGTC PYRGLPAFKP GEKIALPGNF LAADM
Source	Baculovirus
Target Names	H0403D02.15
Protein Names	Recommended name: Anamorsin homolog 1 Alternative name(s): Fe-S cluster assembly protein DRE2 homolog 1
Expression Region	1-265
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.