



Recombinant *Oryza sativa* subsp. *indica* Proteasome subunit alpha type-4-1 (Osl_021120)

Product Code	CSB-MP386100OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2Y9X7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSRRYDSRTT IFSPEGRLYQ VEYAMEAIGN AGSALGVLAA DGVVLVGEKK VTSKLLQTSR SAEKMYKIDS HLACAVAGIM SDANILLNTA RLHAQRYALS YQEPIPVEQL VQSLCDTKQG YTQFGGLRPF GVSFLFAGWD KHHGFQLYMS DPSGNYSGWK AAAGVANSQA AQSMLKQDYR DGMTREEAVA LALKVLSKTM DSTSLTAEKL ELAEVFLQPG TGEVQYQVCS PEAMGKLLAK AGLSQPAPEA
Source	Mammalian cell
Target Names	Osl_021120
Protein Names	Recommended name: Proteasome subunit alpha type-4-1 EC= 3.4.25.1 Alternative name(s): 20S proteasome alpha subunit C 20S proteasome subunit alpha-3
Expression Region	1-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.