



# Recombinant *Oryza sativa* subsp. *indica* Alpha-amylase isozyme 2A (AMYC2)

<b>Product Code</b>	CSB-EP386105OFF-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A2YGY2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	DKILFQGFN WESWRQSGGW YNLLMGKVDD IVAAGVTHVW LPPPSHSVST QGYMPGRLYD LDASRYGTSM ELKSLISALH GKGIQAIADV VINHRCADYK DSRGIYCIFE GGTPDGRLDW GPHMICRDDT QFSDGTGNLD TGADFAAAPD IDHLNGVVQR ELTDWLLWLK SDEVGFDWR LDFARGYSPE VAKVYIEGTT PVGLAVAELW DSMAYGGDGK PEYNQDAHRQ ALVDWVDRVG GTASAGMVFD FTTKGIMNTA VEGELWRLID QQGKAPGVIG WWPAAKAVTFV DNHDTGSTQQ MWPFP SDKVM QGYAYILTHP GNPCIFYDHF FDWGLKEQIA ALVAVRQRNG VTATSSLKIM LHDADAYVAE IDGKVVMKIG SRYDVSSLIP PGFHAAHGN GYAVWEKSAA AAAAAADHRT SSSASL
<b>Source</b>	E.coli
<b>Target Names</b>	AMYC2
<b>Protein Names</b>	Recommended name: Alpha-amylase isozyme 2A EC= 3.2.1.1 Alternative name(s): 1,4-alpha-D-glucan glucanohydrolase Alpha-amylase isozyme C2
<b>Expression Region</b>	22-446
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.